

Updated 01/18/2022 9:17 a.m.

SAMPLE	CODE	NAME	METHOD	ACCREDITED
<b>PRICK TEST EXTERNO</b>				
-	WHITE_FISH_PRICK	PRICK TEST WHITE FISH	(PRICK TEST)	✘
-	APPLE_PRICK	PRICK TEST APPLE	(PRICK TEST)	✘
-	COD_PRICK	PRICK TEST COD	(PRICK TEST)	✘
-	COW_MILK_PRICK	PRICK TEST COW-MILK	(PRICK TEST)	✘
-	EGG_PRICK	PRICK TEST EGG	(PRICK TEST)	✘
-	HAZELNUT_PRICK	PRICK TEST HAZELNUT	(PRICK TEST)	✘
-	MELON_PRICK	PRICK TEST MELON	(PRICK TEST)	✘
-	PEACH_PRICK	PRICK TEST PEACH	(PRICK TEST)	✘
-	PEACH_PRUP3_PRICK	PRICK TEST PEACH PEEL (PRU P 3)	(PRICK TEST)	✘
-	PEANUT_PRICK	PRICK TEST PEANUT	(PRICK TEST)	✘
-	SHRIMP_PRICK	PRICK TEST SHRIMP	(PRICK TEST)	✘
-	WALNUT_PRICK	PRICK TEST WALNUT	(PRICK TEST)	✘
<b>UNICAP EXTERNA</b>				
-	F13_EXTERNO	PEANUT EXTERNAL	(INMUNOCAP)	✘
-	F17_EXTERNO	HAZELNUT EXTERNAL	(INMUNOCAP)	✘
-	F233_EXTERNO	NGAL D 1, OVOMUCOID, EGG EXTERNAL	(INMUNOCAP)	✘
-	F24_EXTERNO	SHRIMP EXTERNAL	(INMUNOCAP)	✘
-	F256_EXTERNO	WALNUT EXTERNAL	(INMUNOCAP)	✘
-	F2_EXTERNO	MILK EXTERNAL	(INMUNOCAP)	✘
-	F416_EXTERNO	RTRI A 19, WHEAT EXTERNAL	(INMUNOCAP)	✘
-	F420_EXTERNO	RPRU P 3, PEACH EXTERNAL	(INMUNOCAP)	✘
-	F422_EXTERNO	RARA H 1, PEANUT EXTERNAL	(INMUNOCAP)	✘
-	F423_EXTERNO	RARA H 2, PEANUT EXTERNAL	(INMUNOCAP)	✘
-	F424_EXTERNO	RARA H3, PEANUT EXTERNAL	(INMUNOCAP)	✘
-	F425_EXTERNO	RCOR A 8, HAZELNUT EXTERNAL	(INMUNOCAP)	✘
-	F428_EXTERNO	RCOR A 1, HAZELNUT EXTERNAL	(INMUNOCAP)	✘
-	F434_EXTERNO	RMAL D 1, APPLE EXTERNAL	(INMUNOCAP)	✘
-	F435_EXTERNO	RMAL D 3, APPLE EXTERNAL	(INMUNOCAP)	✘
-	F439_EXTERNO	RCOR A 14, HAZELNUT EXTERNAL	(INMUNOCAP)	✘
-	F440_EXTERNO	NCOR A 9, HAZELNUT EXTERNAL	(INMUNOCAP)	✘
-	F441_EXTERNO	RJUG R 1, NUT EXTERNAL	(INMUNOCAP)	✘
-	F442_EXTERNO	RJUG R 3, NUT	(INMUNOCAP)	✘
-	F49_EXTERNO	APPLE EXTERNAL	(INMUNOCAP)	✘
-	F77_EXTERNO	NBOS D 5, BETA-LACTOGLOBULIN, MILK EXTERNAL	(INMUNOCAP)	✘
-	F84_EXTERNO	KIWI EXTERNAL	(INMUNOCAP)	✘
-	F95_EXTERNO	PEACH EXTERNAL	(INMUNOCAP)	✘
<b>SÍNTOMAS</b>				